

Premier Products for Women's Healthcare

# Clinic & Practice-Based Products Catalog



Helping you create  
**stronger relationships** with patients



CooperSurgical®



We are committed to supporting healthcare professionals in strengthening the vital connection between women and their OB/Gyn. Our comprehensive portfolio supports the patient journey, from contraception and stem cell planning to diagnosing and treating conditions such as abnormal uterine health, fertility challenges, and pelvic floor disorders, including incontinence.

By partnering with you, we strive to enable personalized care, improve outcomes, and enhance the quality of life for every patient.

## Clinic and Practice-Based Brands

CooperSurgical® is empowering office-based care with an expanded portfolio of innovative devices and products. Within this catalog, you will find an extensive range of instruments and devices that can be used in your practice to improve patient care.

Endosee® Advance

Milex®

LEEP PRECISION™

Fischer Cone Biopsy Excisor®

Pipelle®

Leisegang™

ER-Spec®

Her-Spec®

Euro-Med®

H/S Elliptosphere®

Paragard® Intrauterine  
Copper ContraceptiveCord Blood Registry<sup>SM</sup>

Wallach®

## Table of Contents

### In-Office Hysteroscopy

Endosee® Advance ..... 6

### Milex® Pessaries

Milex® Pessaries for Prolapse and Incontinence ..... 8-11

### LEEP

LEEP PRECISION™ Integrated System ..... 12-13

LEEP Electrodes ..... 14-15

Fischer Cone Biopsy Excisor® ..... 16

LEEP Procedure Essentials ..... 17

LEEP System Accessories ..... 18-19

Euro-Med® LEEP Instruments ..... 20-21

Coated Retractors ..... 22

LEEP Tenacula, Hooks and Rulers ..... 24-25

Wallach® Electrosurgical System &amp; Accessories ..... 26-27

### Colposcopy

Leisegang OptiK® Series ..... 28-31

OptiK Model 1 ..... 28-29

OptiK Model 2 ..... 30-31

Wallach® Zoom Magnification Colposcopes ..... 32

### Gynecological Specialty Instruments

Lighted Specula – Vaginal and Anal ..... 34-35

Euro-Med® Stainless Steel Biopsy Punches ..... 36-37

Euro-Med® Rotating Biopsy Punch Handles ..... 38

Euro-Med® Standard Biopsy Punches ..... 38-39

Euro-Med® Stainless Steel Specula ..... 40-43

Euro-Med® Stainless Steel Endocervical Specula ..... 44-45

Euro-Med® Stainless Steel Lateral Wall Retractors ..... 46

Euro-Med® Stainless Steel Uterine Sounds and Dilators ..... 48-49

Euro-Med® Stainless Steel Curettes ..... 50-51

Euro-Med® Stainless Steel Hooks ..... 52

Euro-Med® Stainless Steel Cervical Tenacula ..... 53

Euro-Med® Stainless Steel Forceps ..... 54-56

Euro-Med® Stainless Steel Needle Holders ..... 56

Euro-Med® Stainless Steel Scissors ..... 57

Endometrial Samplers ..... 58

Reusable Dilators ..... 60

Disposable Sound and Dilators ..... 61

### General Examination: Diagnostic Disposables

Cervical and Endocervical Samplers ..... 62

Diagnostic Disposables ..... 63

### Cervical Cryosurgery

Leisegang™ LM-900 Cryosurgery System ..... 64-65

Cryosurgical Systems ..... 66-67

### Fertility and Family Planning

H/S Elliptosphere Catheters ..... 68

Fertility Assessment ..... 68

H/S Catheters ..... 69

Paragard® Intrauterine Copper Contraceptive ..... 70-71

CBR® Newborn Stem Cell Preservation ..... 72

### Obstetrical Care

Fetal Monitors ..... 73

Obstetrical Dopplers and Probes ..... 74-75

### Vascular Assessment

Vascular Dopplers and Probes ..... 76-77

ABI Vascular Systems ..... 78-81

### Dermatologic Cryosurgery

Ophthalmic Cryosurgery ..... 82

UltraFreeze™ ..... 83

Index/Customer Service ..... 84-85

Warranty Information ..... 85

☑ Available for sale where CE mark is required

🇨🇦 Available for sale in Canada

# Endosee® Advance

## In-Office Diagnostic Hysteroscopy Device

Direct visualization system for diagnostic and therapeutic procedures.



### Endosee® Advance

NUMBER	DESCRIPTION	CA+ CE
ES9000	Endosee Advance Display Module, 1/Box	CA+ CE
ESPX5	Endosee Advance Cannula, Compatible with ES9000 only, 5/Box	CA+ CE
ES-TRAY	Endosee Advance Convenience Kit, 5/Box	
ES-TRAY-IV	Endosee Advance IV Convenience Kit, 5/Box	
<b>REUSABLE INSTRUMENTS (STAINLESS STEEL)</b>		
ES-BPSY	Biopsy Forceps, 5 Fr, 455 mm, 1/Box	
ES-LNGR	Spoon Forceps, Long Serrated, 5 Fr, 455 mm, 1/Box	
ES-SCIS	Scissors, Single Action, 5 Fr, 455 mm, 1/Box	
<b>SINGLE-USE INSTRUMENTS (STERILE)</b>		
ES-DISP-AGR	Disposable Alligator Grasper Forceps, 5 Fr, 460 mm, 5/Box	
ES-DISP-BPY	Disposable Biopsy Forceps, 5 Fr, 460 mm, 5/Box	
ES-DISP-SCS	Disposable Single Action Scissors, 5 Fr, 460 mm, 5/Box	
ES-DISP-10-BL	Disposable Retrieval Snare, 5 Fr, 465 mm, 5/Box	
ES-DISP-15-BL	Disposable Retrieval Snare, 5 Fr, 473 mm, 5/Box	

Biopsy Forceps  
ES-BPSY



Spoon Forceps, Long Serrated  
ES-LNGR



Scissors, Single Action  
ES-SCIS



Disposable Alligator  
Grasper Forceps  
ES-DISP-AGR



Disposable  
Biopsy Forceps  
ES-DISP-BPY



Disposable  
Single Action Scissors  
ES-DISP-SCS



Disposable  
Retrieval Snare  
ES-DISP-10-BL  
ES-DISP-15-BL



To place an order, call us at: 800.243.2974 or visit us online at: [www.coopersurgical.com](http://www.coopersurgical.com)

# Milex® Pessaries for Prolapse & Incontinence

P = Prolapse I = Incontinence



## Parts & Accessories

NUMBER	DESCRIPTION	IN.	MM
MXFIT0012 CA* CE	Fitting Gellhorn Flexible	2 1/4	57
MXFIT0013 CA* CE	Fitting Gellhorn Flexible	2 1/2	64
MXFIT0014 CA* CE	Fitting Gellhorn Flexible	2 3/4	70
MXFIT0015 CA* CE	Fitting Gellhorn Flexible	3	76
MXFIT0016 CA*	Fitting Gellhorn Short Stem Flexible	2 1/4	57
MXFIT0017 CA*	Fitting Gellhorn Short Stem Flexible	2 1/2	64
MXFIT0018 CA*	Fitting Gellhorn Short Stem Flexible	2 3/4	70
MXFIT0019 CA*	Fitting Gellhorn Short Stem Flexible	3	76
MXFIT0020 CA* CE	Fitting Ring with support #2	2 1/4	57
MXFIT0021 CA* CE	Fitting Ring with support #3	2 1/2	64
MXFIT0022 CA* CE	Fitting Ring with support #4	2 3/4	70
MXFIT0023 CA* CE	Fitting Ring with support #5	3	76
MXFIT0024 CA* CE	Fitting Ring with support #6	3 1/4	83
MXFIT0025 CA* CE	Fitting Ring with support #7	3 1/2	89

NUMBER	DESCRIPTION	IN.	MM
MXFIT0801 CA* CE	Fitting Incont. Dish with support #2	2 9/16	65
MXFIT0802 CA* CE	Fitting Incont. Dish with support #3	2 3/4	70
MXFIT0803 CA* CE	Fitting Incont. Dish with support #4	2 15/16	75
MXFIT0804 CA* CE	Fitting Incont. Dish with support #5	3 3/8	80
MXFIT0805 CA* CE	Fitting Donut #2	2 1/2	64
MXFIT0806 CA* CE	Fitting Donut #3	2 3/4	70
MXFIT0807 CA* CE	Fitting Donut #4	3	76
MXFIT0808 CA* CE	Fitting Donut #5	3 3/4	83
MXFIT0809 CA* CE	Fitting Cube #2	1 3/8	35
MXFIT0810 CA* CE	Fitting Cube #3	1 1/2	38
MXFIT0811 CA* CE	Fitting Cube #4	1 5/8	41
MXFIT0812 CA* CE	Fitting Cube #5	1 3/4	44

## Pessary Fitting Kit

NUMBER	DESCRIPTION
MXFIT-CO	Pessary Fitting Set Case and Tray only (no pessaries)
MXFIT-CA CA* CE	Pessary Fitting Set

## A. Ring

For mild prolapse.  
Product # MXPER (Pessary Only)  
**P**

### Ring

NUMBER	IN.	MM	NUMBER	IN.	MM
MXPER00	1 1/4	44	MXPER06	3 1/4	83
MXPER01	2	51	MXPER07	3 1/2	89
MXPER02	2 1/4	57	MXPER08	3 3/4	95
MXPER03	2 1/2	64	MXPER09	4	102
MXPER04	2 3/4	70	MXPER10	4 1/4	108
MXPER05	3	76			

## B. Ring with Knob

For mild prolapse complicated by stress urinary incontinence.  
Product # MXPERK (Pessary Only)  
**P, I**

### Ring with Knob

NUMBER	IN.	MM	NUMBER	IN.	MM
MXPERK00	1 1/4	44	MXPERK06	3 1/4	83
MXPERK01	2	51	MXPERK07	3 1/2	89
MXPERK02	2 1/4	57	MXPERK08	3 3/4	95
MXPERK03	2 1/2	64	MXPERK09	4	102
MXPERK04	2 3/4	70	MXPERK10	4 1/4	108
MXPERK05	3	76			



## C. Ring with Support

For mild prolapse complicated by a mild cystocele.

Product # MXPRS (Pessary Only)  
**P**

### Ring with Support

NUMBER	IN.	MM
MXPRS00	1 1/4	44
MXPRS01	2	51
MXPRS02	2 1/4	57
MXPRS03	2 1/2	64
MXPRS04	2 3/4	70
MXPRS05	3	76
MXPRS06	3 1/4	83
MXPRS07	3 1/2	89
MXPRS08	3 3/4	95
MXPRS09	4	102
MXPRS10	4 1/4	108

## D. Incontinence Ring

For stress urinary incontinence.

Product # MXPCON (Pessary Only)  
**I**

### Incontinence Ring

NUMBER	IN.	MM
MXPCON00	1 1/4	44
MXPCON01	2	51
MXPCON02	2 1/4	57
MXPCON03	2 1/2	64
MXPCON04	2 3/4	70
MXPCON05	3	76
MXPCON06	3 1/4	83
MXPCON07	3 1/2	89
MXPCON08	3 3/4	95
MXPCON09	4	102
MXPCON10	4 1/4	108

## G. Gellhorn

For 3rd degree prolapse/procidentia.  
Product # MXPGE (Pessary Only)

Also available in:  
Short Stem Product # MXPGSS  
Rigid Product # MXPGRS (Pessary Only)  
Short Stem  
**P**

### Gellhorn

NUMBER	IN.	MM
MXPGE1-1/2	1 1/2	38
MXPGE1-3/4	1 3/4	44
MXPGE2-	2	51
MXPGE2-1/4	2 1/4	57
MXPGE2-1/2	2 1/2	64
MXPGE2-3/4	2 3/4	70
MXPGE3-	3	76
MXPGE3-1/4	3 1/4	83
MXPGE3-1/2	3 1/2	89
MXPGE3-3/4	3 3/4	95

## E. Inflatoball™

Inflatable for 3rd degree prolapse/procidentia with or without a mild cystocele and/or rectocele.

Product # MXPINF (Pessary Only)  
**P**

### Inflatoball-Size

NUMBER	IN.	MM
MXPINF5	2	51
MXPINFM	2 1/4	57
MXPINFL	2 1/2	64
MXPINFXL	2 3/4	70

## H. Ring w/ Support & Knob

For mild prolapse complicated by a mild cystocele and stress urinary incontinence.

Product # MXPRSK (Pessary Only)  
**P, I**

### Ring w/ Support & Knob

NUMBER	IN.	MM
MXPRSK00	1 1/4	44
MXPRSK01	2	51
MXPRSK02	2 1/4	57
MXPRSK03	2 1/2	64
MXPRSK04	2 3/4	70
MXPRSK05	3	76
MXPRSK06	3 1/4	83
MXPRSK07	3 1/2	89
MXPRSK08	3 3/4	95
MXPRSK09	4	102
MXPRSK10	4 1/4	108

## F. Shaatz

For mild prolapse complicated by a mild cystocele.

Product # MXPSH (Pessary Only)  
**P**

### Shaatz

NUMBER	IN.	MM
MXPSH1-1/2	1 1/2	38
MXPSH1-3/4	1 3/4	44
MXPSH2-	2	51
MXPSH2-1/4	2 1/4	57
MXPSH2-1/2	2 1/2	64
MXPSH2-3/4	2 3/4	70
MXPSH3-	3	76
MXPSH3-1/4	3 1/4	83
MXPSH3-1/2	3 1/2	89
MXPSH3-3/4	3 3/4	95

## I. Donut

Support for 3rd degree prolapse/procidentia. Not to be deflated during insertion or removal. See Inflatoball if a deflatable pessary is necessary.






Product # MXPDO (Pessary Only)  
**P**

### Donut

NUMBER	IN.	MM
MXPDO0	2	51
MXPDO01	2 1/4	57
MXPDO02	2 1/2	64
MXPDO03	2 3/4	70
MXPDO04	3	76
MXPDO05	3 1/4	83
MXPDO06	3 1/2	89
MXPDO07	3 3/4	95

# Milex® Pessaries for Prolapse & Incontinence

P = Prolapse I = Incontinence



## A. Incontinence Dish

For stress urinary incontinence with mild prolapse.

Product # MXPCOND (Pessary Only)

**P, I**

## B. Incontinence Dish w/Support

For stress urinary incontinence with mild prolapse and a mild cystocele.

Product # MXPCONDS (Pessary Only)

**P, I**

### Incontinence Dish

Dish		Dish w/ Support	
NUMBER	NUMBER	IN.	MM
MXPCOND00	MXPCONDS00	2 1/8	55
MXPCOND01	MXPCONDS01	2 3/8	60
MXPCOND02	MXPCONDS02	2 1/2	65
MXPCOND03	MXPCONDS03	2 3/4	70
MXPCOND04	MXPCONDS04	3	75
MXPCOND05	MXPCONDS05	3 1/8	80
MXPCOND06	MXPCONDS06	3 3/8	85

## C. Gehrung with Knob/Folding

For cystocele with or without mild rectocele complicated by stress urinary incontinence.

Product # MXPGSK (Pessary Only)

**I**

## D. Gehrung/Folding

Support for cystocele and rectocele. It is also more adaptable in cases of procidentia where the uterus tends to herniate when other pessaries are used.

Product # MXPGS (Pessary Only)

**P**

### Gehrung

Gehrung w/ Knob/Folding		Gehrung Folding					
NUMBER	NUMBER	IN. (A)	MM (A)	IN. (B)	MM (B)	IN. (C)	MM (C)
MXPGSK00	MXPGS00	1 1/8	27	1 1/2	38	1 3/8	35
MXPGSK01	MXPGS01	1 1/8	29	1 3/4	44	1 1/2	38
MXPGSK02	MXPGS02	1 3/8	35	2	51	1 5/8	41
MXPGSK03	MXPGS03	1 1/2	38	2 1/8	54	1 3/4	44
MXPGSK04	MXPGS04	1 5/8	41	2 1/4	57	1 7/8	48
MXPGSK05	MXPGS05	1 3/4	44	2 3/8	60	2	51
MXPGSK06	MXPGS06	1 7/8	48	2 1/2	64	2 1/8	54
MXPGSK07	MXPGS07	2	51	2 7/8	73	2 1/4	57
MXPGSK08	MXPGS08	2 1/8	54	3	76	2 3/8	60

## E. Cube

Support for 3rd degree prolapse/procidentia. Should be removed nightly.

Product # MXPEC (Pessary Only)

Also available in:  
Cube with Drainage Holes  
Product # MXPECH (Pessary Only)

**P, I**

### Cube

Cube		w/ Holes	
NUMBER	NUMBER	IN.	MM
MXPEC00	MXPECH00	1	25
MXPEC01	MXPECH01	1 3/16	30
MXPEC02	MXPECH02	1 3/8	35
MXPEC03	MXPECH03	1 1/2	38
MXPEC04	MXPECH04	1 5/8	41
MXPEC05	MXPECH05	1 3/4	44
MXPEC06	MXPECH06	2	51
MXPEC07	MXPECH07	2 1/4	57

# LEEP PRECISION™ Integrated System

Loop Electrosurgical Excision Procedures



Fully integrated and mobile, the LEEP PRECISION System combines three components – Cart, Generator and Smoke Evacuator.



## LEEP PRECISION Integrated System

NUMBER	DESCRIPTION	
LP-10-120	LEEP PRECISION Integrated System, 120V	CA*
LP-10-220	LEEP PRECISION Integrated System, 220V International	CE

## LEEP PRECISION Generator

NUMBER	DESCRIPTION	
LP-20-120	LEEP PRECISION Generator, 120V	CA*
LP-20-220	LEEP PRECISION Generator, 220V International	CE

## LEEP PRECISION Smoke Evacuator

NUMBER	DESCRIPTION	
LP-30-120	LEEP PRECISION Smoke Evacuator, 120V, Stand-Alone	CA*
LP-30-220	LEEP PRECISION Smoke Evacuator, 220V, Stand-Alone International	CE

## LEEP PRECISION Cart

NUMBER	DESCRIPTION	
LP-40-001	LEEP PRECISION Cart	CA* CE

## LEEP PRECISION Generator Disposables

NUMBER	DESCRIPTION	
LP-50-101	Hand Switch-Operated Push-Button Cautery Pencils (10/Box)	CA*
LP-50-102	Foot Pedal-Operated Cautery Pencils (10/Box)	CA*
LP-50-103	Foot Pedal Adapter (1/Bag)	CA*
LP-50-201	Dual Dispersive Electrodes with Connector (10/Box)	CA*
LP-50-203	Reusable Cable (1/Bag)	CA*



## LEEP PRECISION Smoke Evacuator Disposables

NUMBER	DESCRIPTION	
52560	Smoke Evacuator Pre-Filter, Single-Patient Use (10/Bag)	CA*
6082	Smoke Evacuator ULPA Filter (1/Box)	CA*
6083	Smoke Evacuator Reducer, Single-Patient Use (10/Bag)	CA*
6084	Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)	CA*

# LEEP Electrodes

(Shown Actual Size)

Our broad selection of 21 different offerings includes loops, squares, needle and ball tips specifically designed for loop electrosurgical excision procedures.



## Large-Radius Loop Electrodes ≥ 2 cm Wide CA\* CE

NUMBER	DESCRIPTION
R2510	2.5 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R2020	2.0 cm W x 2.0 cm D • 12 cm shaft (5/Box)
R2015	2.0 cm W x 1.5 cm D • 12 cm shaft (5/Box)
R2010	2.0 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R2008	2.0 cm W x 0.8 cm D • 12 cm shaft (5/Box)

## Medium-Radius Loop Electrodes <2 cm Wide CA\* CE

NUMBER	DESCRIPTION
R1515	1.5 cm W x 1.5 cm D • 12 cm shaft (5/Box)
R1507	1.5 cm W x 0.7 cm D • 12 cm shaft (5/Box)
R1505	1.5 cm W x 0.5 cm D • 12 cm shaft (5/Box)
R1010	1.0 cm W x 1.0 cm D • 12 cm shaft (5/Box)
R1007	1.0 cm W x 0.7 cm D • 12 cm shaft (5/Box)

## Square Electrodes CA\* CE

NUMBER	DESCRIPTION
S1010	1.0 cm W x 1.0 cm D 12 cm shaft (5/Box)
S1007	1.0 cm W x 0.7 cm D 5.5 cm shaft 15° bend (5/Box)
S1004	1.0 cm W x 0.4 cm D 5.5 cm shaft 15° bend (5/Box)
V1004	1.0 cm W x 0.4 cm D 8 cm shaft (5/Box)

## Needle Electrodes CA\* CE

NUMBER	DESCRIPTION
N1200	0.8 mm D x 16 mm 12 cm shaft (5/Box)
N0550	0.8 mm D x 16 mm 5.5 cm shaft (5/Box)
N2512	0.25 mm D x 8 mm 12 cm shaft (5/Box)
N2555	0.25 mm D x 5 mm 5.5 cm shaft (5/Box)

## Ball Electrodes CA\* CE

NUMBER	DESCRIPTION
B0512	5 mm D 12 cm shaft (5/Box)
B0555	5 mm D 5.5 cm shaft (5/Box)
B0312	3 mm D 12 cm shaft (5/Box)


# Fischer Cone Biopsy Excisor<sup>®</sup>

For Cone Biopsy Procedures

The Fischer Cone Biopsy Excisor (FCBE) comes in a variety of sizes and shapes.

## Fischer Cone Biopsy Excisor CA\* CE

NUMBER	DESCRIPTION
900-150	Fischer Cone Biopsy Excisor, Small (5/Box)
900-151	Fischer Cone Biopsy Excisor, Medium (5/Box)
900-152	Fischer Cone Biopsy Excisor, Large (5/Box)
900-154	Fischer Cone Biopsy Excisor, Medium Extended (5/Box)
900-155	Fischer Cone Biopsy Excisor, Large Shallow (5/Box)
900-156	Fischer Cone Biopsy Excisor, Mixed Set (S, M, ME, L & LS) 1 each (5/Box)
900-157	Fischer Cone Biopsy Excisor, Small Wide-Angle (5/Box)
900-158	Fischer Cone Biopsy Excisor, Small Wide-Angle Extended (5/Box)



# LEEP Procedure Essentials

## Potocky Needle<sup>®</sup>

27-gauge Potocky Needle. The Potocky Needle attaches to a "dental style" cartridge syringe, while the Potocky Needle with Luer Lock is used with a standard medical syringe.

NUMBER	DESCRIPTION
6066	Potocky Needles (6/Box)
6166	Potocky Needles with Luer Lock (6/Box)



## Needle Extension

NUMBER	DESCRIPTION
433-132	Needle Extension • 10 cm (1/Box)
433-133	Needle Extension • 13 cm (1/Box)
433-134	Needle Extension • 15 cm (1/Box)



# LEEP System Accessories



### Foot Pedal, LEEP System 1000

Featuring pneumatic, air-operated control, the foot pedal does not use electric current.

NUMBER	DESCRIPTION
6032	Foot Pedal for LEEP System 1000




### Return Pad, Single-Patient Use CA\*

Non-sterile patient return pads.

NUMBER	DESCRIPTION
6050P1	LEEP Return Pad, Single-Patient Use (10/Box)


# LEEP Accessories



### Compatible Disposable Accessories

See your local sales representative for more information.

NUMBER	DESCRIPTION
395-210	Electrode Tip Cleanser, Disposable (24/Box)
395-302	Conductive Electrode Pads, Disposable for ESU (25/Box)
395-305	Conductive Electrode Pads, Dual Dispersive, Disposable (25/Box) Plus II Gen.
395-310	Hand-Switch Pencil and Cord, Disposable (25/Box)
72308	Dual Section Grounding Pads (10/Box)
72102	Filter/Hose Set for Smoke Evacuator (Pack of 3)
72311	Patient Return Cable



### Tubing Set, Single-Patient Use CA\*

Each set includes 6 feet of sterile disposable tubing and a DSE (Disposable Smoke Evacuator) tube. This set completes the connection between the reducer and the specula.


NUMBER	DESCRIPTION
6084	Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box)



### Handpiece, Single-Patient Use CA\*

The sterile, single-patient use handpiece works with any of our disposable electrodes. The wand is operated with a footswitch.

NUMBER	DESCRIPTION
6040	Handpiece, Single-Patient Use (10/Box)



### Disposable Smoke Evacuator (DSE) Tubing

Non-sterile single-use, DSE tubing is designed to work specifically with the LEEP rose-coated specula. The 9-inch DSE tubing is flexible and noncombustible and contains no harmful PVC.

Also available as a set with 6 feet of tubing. See Tubing Set, Single-Patient Use, product item number 6084.

NUMBER	DESCRIPTION
6150	Disposable Smoke Evacuator Tubing – 9" straw (50/Bag)

### Patient Return Pad Adapter LEEP System 1000



This reusable replacement adapter connects the single-use Patient Return Pad to the LEEP System 1000.

NUMBER	DESCRIPTION
23635	Patient Return Pad Adapter for LEEP System 1000 (for use with 6050P1 pads)



### Reducer CA\*

Disposable, single-patient use reducer connects the Pre-Filter (52560) with the Patient Tubing Set (6084) on the Smoke Evacuator.

NUMBER	DESCRIPTION
6083	Smoke Evacuator Reducer (10/Box)



### Pre-Filter CA\*

Replacement filters for the CooperSurgical Smoke Evacuator. Used to eliminate smoke odor.

NUMBER	DESCRIPTION
52560	Smoke Evacuator Pre-Filter (10/Box)



### ULPA Filter CA\*

Used with the Smoke Evacuator.

NUMBER	DESCRIPTION
6082	Smoke Evacuator ULPA Filter (1/Box)

# Euro-Med® LEEP Instruments

These nonconductive, fully autoclavable LEEP instruments are offered in a variety of styles to facilitate cervical tissue excision and removal.



**Pederson LEEP Specula**

NUMBER	DESCRIPTION	SIZE MM
F223	Small	W13 L81
F221	Medium	W22 L114
F220	Large	W25 L120



**Vu-More™ LEEP Specula**

NUMBER	DESCRIPTION	SIZE MM
F217	Medium	W33 L113
F215	Large	W35 L127



**Vu-Max™ LEEP Specula**

NUMBER	DESCRIPTION	SIZE MM
F222	Pederson Medium	W22 L114
F224	Pederson Large	W25 L120
F218	Graves Medium	W33 L113
F219	Graves Large	W35 L127



**Graves LEEP Specula**

NUMBER	DESCRIPTION	SIZE MM
F205	Medium	W33 L113
F200	Large	W35 L127



**Snowman Specula™**

NUMBER	DESCRIPTION	SIZE MM
F107	Graves Blade	W35 L152
F117	Pederson Blade	W25 L152



**Mini Vu-More™ LEEP Specula**

NUMBER	SIZE MM
F250	W34 L93



**Open-Sided LEEP Specula**

NUMBER	SIZE MM
F210	W35 L107



**Auvard Surgical Specula**

NUMBER	DESCRIPTION	SIZE MM
F230	Medium (Wide/ Short Blade)	W43 L55
F232	Large (Narrow/Long Blade)	W35 L90




**Pederson LEEP Specula (with Integrated Smoke Tube)**

NUMBER	DESCRIPTION	SIZE MM
61521	Small	W16 L84
61518	Medium	W22 L114
61519	Large	W26 L121




**Graves LEEP Specula (with Integrated Smoke Tube)**

NUMBER	DESCRIPTION	SIZE MM
61513	Small	W22 L87
61514	Medium	W33 L113
61515	Large	W35 L125



(W=Distal End) (L=Penetration Length)

# Euro-Med® LEEP Instruments



**Tru-View™ LEEP  
Lateral Wall Retractor**

NUMBER	DESCRIPTION	SIZE MM
F317	Lateral Wall Retractor	W32 L190



**Jackson LEEP Retractor**

NUMBER	DESCRIPTION	SIZE MM
F235	Medium	W40 L89
F234	Large	W40 L108



**Cer-View™ LEEP  
Lateral Wall Retractor**

NUMBER	DESCRIPTION	SIZE MM
F410	Medium	W22 L142
F420	Large	W22 L155




# Euro-Med® LEEP Tenacula, Hooks and Rulers



**LEEP Ring Forceps**

NUMBER	INSTRUMENT LENGTH
F920	220 mm




**LEEP Tissue Forceps**

NUMBER	INSTRUMENT LENGTH
F700	210 mm



**LEEP Campion Forceps**

NUMBER	INSTRUMENT LENGTH
F925	220 mm



**LEEP Single-Tooth Tenaculum**

NUMBER	INSTRUMENT LENGTH
F900	240 mm



**LEEP Emmett Tenaculum**

NUMBER	INSTRUMENT LENGTH
F805	250 mm



**LEEP Two-Prong Hook**

NUMBER	INSTRUMENT LENGTH
F802	260 mm



**LEEP Three-Prong Hook**

NUMBER	INSTRUMENT LENGTH
F803	260 mm



**LEEP Iris Hook**

NUMBER	INSTRUMENT LENGTH
F800	263 mm



**LEEP Hook**

NUMBER	INSTRUMENT LENGTH	
F810	Straight	250 mm
F812	Left-Angled	250 mm
F813	Right-Angled	250 mm



**LEEP Depth Ruler**

NUMBER	INSTRUMENT LENGTH
F600	240 mm



**LEEP Width Ruler**

NUMBER	INSTRUMENT LENGTH
F610	210 mm

## LEEP Endocervical Specula



**Large**


NUMBER	TIP SIZE	LENGTH
F500	5 mm	273 mm



**Medium**

NUMBER	TIP SIZE	LENGTH
F510	4 mm	273 mm

# Electrosurgical System Accessories for Wallach® Quantum 2000™



### Foot Pedal–Operated Pencils

Sterile, single-use only. For use with Quantum 2000™ and other comparable electrosurgical generators. Reusable adapter required. Use with any standard disposable electrodes.


NUMBER	DESCRIPTION
909028	Foot Pedal–Operated Pencils, Sterile (25 per Box)



### Disposable Smoke Evacuator Tubing with Attached Speculum Tubing

Sterile, single-use six-foot-long tubing connects to the smoke evacuator filter assembly.

NUMBER	DESCRIPTION
920002	Disposable Smoke Evacuator Tubing with Speculum Tubing (10 per Box)




### Disposable Smoke Evacuator Speculum Tubing (Non-Sterile)

NUMBER	DESCRIPTION
920001	Disposable Speculum Tubing (50 per Box)

### Disposable Split Grounding Pads (with Connector)

For use with Quantum 2000™.

NUMBER	DESCRIPTION
909078	Disposable Split Grounding Pads with Connector (10 per Box)



# Leisegang OptiK® Model 1

Optical Devices Designed for Non Contact Visualization of the Tissues of the Cervix, Vulva, Vagina, And Anogenital Areas



High-Intensity LED Illumination



## LEISEGANG OptiK® Model 1

NUMBER	DESCRIPTION
OPTIK1-00	Leisegang OptiK Model 1 Tilt Colposcope (1DL LED): 3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, Tilt Base
OPTIK1-01	Same as OPTIK1-00 except with Upright Rolling Base
OPTIK1-02	Same as OPTIK1-00 except with medium-weight Swing-Arm Rolling Base
OPTIK1-03	Same as OPTIK1-00 except with heavy-weight Balance-O-Matic Rolling Base

## Four Base Choices:



Tilt Base  
OPTIK1-00



Upright Rolling Base  
OPTIK1-01



Swing-Arm Rolling Base  
OPTIK1-02



Balance-O-Matic Rolling Base  
OPTIK1-03

# Leisegang OptiK® Model 2


OptiK Model 2 with  
a Digital SLR Camera



OptiK Model 2 with  
an Integrated Camera

## Three Base Choices:


OptiK Model 2 is available in three base configurations including Upright, Swing-Arm and Balance-O-Matic. The OptiK Model 2 features a working distance of 300 mm, field of view of 19 mm to 76 mm, and depth of field of 2 mm to 8 mm.



Upright Rolling Base  
OPTIK2-03-002



Swing-Arm Rolling Base  
OPTIK2-04-002



Balance-O-Matic Rolling Base  
OPTIK2-05-002

## Leisegang OptiK® Model 2 Sales Configuration Options:

### LEISEGANG OptiK® Model 2, Base and dSLR Camera CA\* CE

NUMBER	DESCRIPTION
OPTIK2-03-002	3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, 1.0" CCD optimized, medium-weight Upright Rolling Base, T2 adapter/tubus, Canon camera kit, battery, charger, USB wire, A/C adapter, wired remote and memory card
OPTIK2-04-002	Same as OPTIK2-03-002 except with Swing Arm Rolling Base
OPTIK2-05-002	Same as OPTIK2-03-002 except with heavy-weight Balance-O-Matic Rolling Base


## OptiK Model 2 Accessories

### Cameras and Monitors

NUMBER	DESCRIPTION
53576	Canon® Camera Kit for OptiK Model 2 1.0" CCD only
23676FG	Foot Pedal for Canon Camera

# Wallach® Zoom Magnification Colposcopes

Continuously Variable Nikon® 4.5 – 20X Zoom, Trulight LED Lighting



## ZoomScope® with Trulight™ CA\*

NUMBER	DESCRIPTION
906043-SP-3	ZoomScope Overhead arm colposcope Quantum series. Smart purchase package: includes a set of sterilizable handles, a set of sterilizable knob covers for zoom and focus, dust cover; 3-leg base and seven year warranty
906043-SP-4	ZoomScope with Trulight and 4-Leg Base
906043-40TU-5	ZoomScope Trulight Video, Digital USB Camera
906057T	ZoomStar with Trulight
906057-40TU	ZoomStar with Trulight Video



ZoomScope® with Trulight™, Overhead Suspension Arm

ZoomStar® with Trulight™, on a Compact Mobile Stand



# Lighted Specula – Vaginal and Anal



## ER-Spec®

Single-use vaginal speculum with integrated LED light source.

NUMBER	DESCRIPTION
C020130-1	ER-Spec® (Extra-Small)
C020100-1	ER-Spec® (Small)
C020110-1	ER-Spec® (Medium)
C020120-1	ER-Spec® (Large)
C020140-1	ER-Spec® (Extra-Large)
C020300	ER-Spec® (Small -Sterile)
C020310	ER-Spec® (Medium - Sterile)
C020320	ER-Spec® (Large - Sterile)



## Her-Spec®

Single-use click-free vaginal speculum with integrated LED light source.

NUMBER	DESCRIPTION
C020500	Her-Spec® Small, Non-Sterile (1 Box of 16)
C020510	Her-Spec® Medium, Non-Sterile (1 Box of 16)
C020520	Her-Spec® Large, Non-Sterile (1 Box of 14)
C020600	Her-Spec® Small, Sterile (1 Box of 10)
C020610	Her-Spec® Medium, Sterile (1 Box of 8)
C020620	Her-Spec® Large, Sterile (1 Box of 8)

## Accessories

NUMBER	DESCRIPTION
C040100	OfficePack® for Hysteroscopy (1 Box of 5)
C050100	Fluid Management Pack (1 Box of 5)
C070100	Vaginal Exam Pack (1 Box of 20)
C070110	IUD Pack (1 Box of 10)



## OfficeSpec®

Single-use right side-opening vaginal speculum with integrated LED light source.

NUMBER	DESCRIPTION
C020001	OfficeSpec® (Right Side-Opening Medium)



## Anospec®

Single-use speculum designed to aid in the examination of the anus to screen, diagnose, evaluate, and treat problems of the anal canal.

NUMBER	DESCRIPTION
C060100	Anospec®, 90 mm x 18 mm slotted, Non-Sterile (1 Box of 20)
C060110	Anospec®, 96 mm x 23 mm slotted, Non-Sterile (1 Box of 16)
C060120	Anospec®, 103 mm x 18 mm beveled, Non-Sterile (1 Box of 20)

# Euro-Med™ Stainless Steel Biopsy Punches

Blue markings indicate bite shape.

## A. Tischler-Morgan™ Biopsy Punch

The standard bite incorporates a single tooth on the upper and lower jaws that draws tissue in for sample retention.

NUMBER	LENGTH	BITE
64-689	21 cm	7 mm x 3 mm
64-689L	25 cm	7 mm x 3 mm

## B. Baby Tischler Biopsy Punch

Small bite incorporates a tooth permitting sampling at the endocervical canal.

NUMBER	LENGTH	BITE
64-691	21 cm	4.1 mm x 2.3 mm
64-691L	26 cm	4.1 mm x 2.3 mm

## C. Tischler-Kevorkian™ Biopsy Punch

Large rectangular bite includes upper and lower teeth.

NUMBER	LENGTH	BITE
64-679	22 cm	9 mm x 3 mm

## D. Kevorkian-Pacific™ Biopsy Punch

A slightly smaller bite and additional anchoring teeth on the bottom jaw.

NUMBER	LENGTH	BITE
64-675	22 cm	9 mm x 3 mm

## E. Mini-Townsend™ Biopsy Punch

NUMBER	LENGTH	BITE
64-690	21 cm	4 mm x 2.3 mm

## F. Coppleson Biopsy Punch



The unique jaw shape yields a distinct sample orientation for histologic analysis.

NUMBER	LENGTH	BITE
64-660	22 cm	7.5 mm x 3 mm

## G. Baggish Biopsy Punch

This punch provides an oval-shaped bite. The jaws have two bottom teeth and a serrated bed.

NUMBER	LENGTH	BITE
64-695	22 cm	6.5 mm x 3.5 mm



## H. Burke™ Biopsy Punch

The design features an open handle and an oval-shaped bite for small samples.

NUMBER	LENGTH	BITE
64-677	22 cm	5.5 mm x 3 mm



## I. Euro-Med™ Keyes Vulvar Punch Set and Singles


This set includes six stainless steel interchangeable tips ranging from 3 mm to 8 mm in diameter, and a universal handle. The individual tips can be placed on the threaded handle as required.

NUMBER	LENGTH
68-180	Set includes: five tips, storage tray and wooden box

# Euro-Med™ Rotating Biopsy Punch Handles

A 360-Degree Rotating Shaft for Sample Collection

These rotating punch tips fit either the Euro-Med™ closed rotating handle or the Townsend open rotating handle.



### Tischler-Morgan™ Tip

NUMBER	LENGTH	BITE
64-484	23 cm	7 mm x 3 mm

### Mini Tip Down

NUMBER	LENGTH	BITE
64-656	23 cm	7 mm x 3 mm

### Tischler Tip Up

NUMBER	LENGTH	BITE
64-751	23 cm	7 mm x 3 mm

### Mini-Townsend™

NUMBER	LENGTH	BITE
64-460	23 cm	4 mm x 2.3 mm

### Tischler Tip Down

NUMBER	LENGTH	BITE
64-756	23 cm	7 mm x 3 mm

### Oval Tip

NUMBER	LENGTH	BITE
64-489	23 cm	5 mm x 3 mm

### Mini Tip Up

NUMBER	LENGTH	BITE
64-651	23 cm	7 mm x 3 mm

### Kevorkian Tip

NUMBER	LENGTH	BITE
64-488	23 cm	8 mm x 3 mm

# Euro-Med™ Standard Biopsy Punches



Handle for Eppendorfer, Wittner and Kevorkian Tips

Handle for Schubert Tip

### Eppendorfer Tip

This oval cutting jaw obtains an 8 mm x 4.5 mm cervical specimen. A basket in the jaw helps retain the biopsy sample.

NUMBER	LENGTH	BITE
64-687	21 cm	8 mm x 4.5 mm

### Wittner Tip

The angle of the jaw enables extra access to areas of the cervix. The large teeth and strong bite provide a maximum tissue sample area.

NUMBER	LENGTH	BITE
64-681	21.5 cm	8 mm x 4 mm

### Kevorkian Tip

This punch facilitates simple insertion and sampling of the cervical canal.

NUMBER	LENGTH	BITE
64-685	26.5 cm	9 mm x 3 mm

### Schubert Tip

Features a large bite and an angled handle.


NUMBER	LENGTH	BITE
64-670	21.5 cm	8 mm x 4.5 mm

# Euro-Med™ Stainless Steel Specula

## A. Graves €

Available with standard yoke and an open-sided configuration. The Graves blade is also available in longer lengths with a Vu-More™ or Vu-Max™ yoke.

NUMBER	DESCRIPTION	SIZE MM
64-100	Small	W19 L82
64-102	Medium	W33 L114
64-104	Large	W35 L127
64-108	Open-sided	W35 L100




A.

## B. Pederson

Pederson vaginal speculum is a double-bladed instrument with a screw mechanism to keep the blades open once inserted into the vaginal canal. The speculum is designed to aid in the visualization of the vaginal canal for diagnostic and/or procedural purposes.

NUMBER	DESCRIPTION	SIZE MM
391-671	Virgin	W13 L78
64-110	Small	W13 L82 €
391-674	Extra-narrow	W17 L108 €
64-112	Medium	W22 L114 €
64-114	Large	W25 L130 €




B.

## C. Vu-More™ €

The Vu-More provides a 25 percent larger viewing area at the yoke.

NUMBER	DESCRIPTION	SIZE MM
64-130	Mini	W34 L93
63-902	Medium	W33 L113
63-900	Large	W35 L127
63-905	X Large	W33 L150
63-906	2X Large	W33 L160
63-907	3X Large	W33 L170
63-908	4X Large	W33 L180

New Longer Length



C.



D.

## D. Vu-Max™ €

Offered in both Pederson and Graves configurations.

NUMBER	DESCRIPTION	SIZE MM
64-113	Pederson Medium	W22 L114
64-115	Pederson Large	W25 L120
64-218	Graves Medium	W33 L114
64-219	Graves Large	W35 L125
64-225	X Large	W35 L150
64-226	2X Large	W35 L160
64-227	3X Large	W35 L170
64-228	4X Large	W35 L180

New Longer Length



E.

## E. Snowman Specula™

Longer blades expose access to the cervix.

NUMBER	DESCRIPTION	SIZE MM
64-107	Graves Blade	W35 L152 €
64-117	Pederson Blade	W25 L152 €
64-127	One each/Graves & Pederson	



G.

## G. Auvard €

This 2.5-pound weighted speculum for posterior vaginal wall retraction in vaginal procedures.

NUMBER	DESCRIPTION	SIZE MM
64-422	Med (Wide/Short Blade)	W43 L55
64-421	Lg (Narrow/Long Blade)	W35 L90




H.

## H. Weisman Graves Side Opening €

Available in two sizes with left or right openings.

NUMBER	DESCRIPTION	SIZE MM
391-667A	Left side open, Med	W33 L102
391-667	Left side open, Lg	W36 L121
391-668A	Right side open, Med	W35 L102
391-668	Right side open, Lg	W36 L120



F.

## F. ExpandaView™ Integrated Speculum and Lateral Retractor

Combines a traditional speculum with a lateral wall retractor.

NUMBER	DESCRIPTION	SIZE MM
64-161	Medium	W22 L109 €
64-162	Large	W22 L133 €
64-163	Stainless Steel Instrument Holder	



I.

## I. Pederson Side Opening €

The Pederson open-sided specula features either a left or right side opening.

NUMBER	DESCRIPTION	SIZE MM
64-112L	Left side open, Med	W22 L114
64-114L	Left side open, Lg	W25 L120
64-112R	Right side open, Med	W22 L114
64-114R	Right side open, Lg	W25 L120

(W=Distal End) (L=Penetration Length)

# Euro-Med™ Stainless Steel Specula



## Specula Blades and Yoke Dimensions

Use our speculum sizing guide to select the proper tools for your practice. Specula are available in a variety of sizes.

See pages 44-45 for a full product selection and model numbers.

### Vu-Max™ Yoke – Graves Blades

Medium – W: 33 mm L: 114 mm  
Large – W: 35 mm L: 125 mm

### Vu-Max Yoke – Pederson Blades

Medium – W: 22 mm L: 114 mm  
Large – W: 25 mm L: 120 mm

### C. Vu-More™ Yoke – Graves Blades

Mini – W: 34 mm L: 93 mm  
Medium – W: 33 mm L: 113 mm  
Large – W: 35 mm L: 127 mm

### D. Standard Yoke – Pederson Blades

Virgin – W: 13 mm L: 78 mm  
Small – W: 13 mm L: 82 mm  
Extra Narrow – W: 17 mm L: 108 mm  
Medium – W: 22 mm L: 114 mm  
Large – W: 25 mm L: 130 mm

Also Available with Side Opening Yoke

### E. Standard Yoke – Graves Blades

Small – W: 19 mm L: 82 mm  
Medium – W: 33 mm L: 114 mm  
Large – W: 35 mm L: 127 mm

Also Available with Side Opening Yoke

### F. Vu-More Yoke – Graves or Pederson Blade

Snowman™ – W: 35 mm L: 152 mm  
Snowman – W: 25 mm L: 152 mm

### G. Vu-Max and Vu-More Yokes – Graves Blades

X Large – W: 35 mm L: 150 mm  
2X Large – W: 35 mm L: 160 mm  
3X Large – W: 35 mm L: 170 mm  
4X Large – W: 35 mm L: 180 mm

# Euro-Med™ Stainless Steel Endocervical Specula



**Ratchet Handle Speculum** CE

Stainless steel endocervical speculum with quick-release lock.

NUMBER	TIP SIZE	LENGTH
66-437	Small 3 mm	10.6"/270 mm
66-435	Medium 4 mm	10.6"/270 mm
66-436	Large 5 mm	10.6"/270 mm



**Screw/Knob Handle Speculum** CE

The screw/knob handle stainless steel speculum adjusts and maintains position using the threaded thumb screw.

NUMBER	TIP SIZE	LENGTH
66-431	Small 3 mm	10.6"/270 mm
66-430	Medium 4 mm	10.6"/270 mm
66-432	Large 5 mm	10.6"/270 mm







**No-Lock Handle** CE

This stainless steel endocervical speculum offers singlehanded action. The No-Lock Handle allows for adjustment as necessary.

NUMBER	TIP SIZE	LENGTH
66-433	Small 3 mm	10.6"/270 mm
67-435	Medium 4 mm	10.6"/270 mm
66-434	Large 5 mm	10.6"/270 mm

Each Endocervical Speculum Available in Small, Medium or Large Tip



# Euro-Med™ Stainless Steel Lateral Wall Retractors



### Cer-View™ Lateral Wall Retractor CA\* CE

NUMBER	DESCRIPTION	SIZE MM
64-310	Medium	W22 L150
64-320	Large	W22 L160

### Tru-View™ Lateral Wall Retractor CE

NUMBER	DESCRIPTION	SIZE MM
64-307	Extra Large	W32 L190

### Jackson Vaginal Retractor CE

NUMBER	DESCRIPTION	SIZE MM
64-322	Medium	W40 L98
64-324	Large	W40 L115



# Euro-Med™ Stainless Steel Uterine Sounds and Dilators

## Dil-Os™ Uterine Sound/Dilator

NUMBER	DESCRIPTION	SIZE	CE
64-336	Slight Curve	L12.7"/322 mm	CE
64-337	Strong Curve	L12.5"/317 mm	CE
64-338	Set of two (64-336, 64-337)		



## Sims Uterine Sound CE

NUMBER	DESCRIPTION	SIZE
64-601	Sims Uterine Sound – Stainless Steel	L12.5"/317 mm
64-600	Sims Uterine Sound – Silver-Plated Copper*	L12.5"/317 mm

\* Silver-plated copper may tarnish.



## Pratt Dilator Set CE

Available individually or in a set of 8.

NUMBER	SIZE	NUMBER	SIZE
391-830	13 Fr/15 Fr – 41 Fr/43 Fr (8/Set)	391-838	33 Fr/35 Fr
391-833	13 Fr/15 Fr	391-839	37 Fr/39 Fr
391-834	17 Fr/19 Fr	391-840	41 Fr/43 Fr
391-835	21 Fr/23 Fr	391-841	45 Fr/47 Fr
391-836	25 Fr/27 Fr	391-842	49 Fr/51 Fr
391-837	29 Fr/31 Fr	391-843	53 Fr/55 Fr
		391-844	57 Fr/59 Fr



## Hegar Dilator Set CE

Stainless steel Hegar dilators provide the option of double-ended dilators. Available individually or in a set of 8.

NUMBER	SIZE MM
61207	9 Fr (3 mm)/12 Fr (4 mm); 51 Fr (17 mm)/ 54 Fr (18 mm) (8/Set)
61208	9 Fr (3 mm)/12 Fr (4 mm)
61209	15 Fr (5 mm)/18 Fr (6 mm)
61210	21 Fr (7 mm)/24 Fr (8 mm)
61211	27 Fr (9 mm)/30 Fr (10 mm)
61212	33 Fr (11 mm)/36 Fr (12 mm)
61213	39 Fr (13 mm)/42 Fr (14 mm)
61214	45 Fr (15 mm)/28 Fr (16 mm)
61215	51 Fr (17 mm)/54 Fr (18 mm)



## Hank Cervical Dilators CE

Hank cervical dilators available individually or in a set of six ranging from 9 French (3 mm) to 20 French (6.66 mm). The set comes in a custom-fitted and numbered pouch.

NUMBER	SIZE MM
64-370	6 Dilators and Storage Pouch
61201	9 Fr (3 mm)/10 Fr (3.33 mm)
61202	11 Fr (3.67 mm)/12 Fr (4 mm)
61203	13 Fr (4.33 mm)/14 Fr (4.67 mm)
61204	15 Fr (5 mm)/16 Fr (5.33 mm)
61205	17 Fr (5.67 mm)/18 Fr (6 mm)
61206	19 Fr (6.33 mm)/20 Fr (6.67 mm)



## Mini Cervical Dilator Set CE

NUMBER	DESCRIPTION
64-374	5 Dilators, Handle and Storage Pouch




## Small Cervical Dilator Set CE

NUMBER	DESCRIPTION
64-372	6 Dilators and Storage Pouch

# Euro-Med™ Stainless Steel Endocervical Curettes

# Euro-Med™ Endometrial Curettes



**Kevorkian with Basket**

NUMBER	LENGTH
66-442	12"/305 mm

**Kevorkian without Basket**

NUMBER	LENGTH
66-440	12"/305 mm

**Townsend™** CE

NUMBER	LENGTH
66-445	12.3"/312 mm

**A. Novak** CE

NUMBER	SIZE	LENGTH
64-624	1 mm	9.8"/248 mm
64-626	2 mm	9.8"/248 mm
64-628	3 mm	9.8"/248 mm
64-630	4 mm	9.8"/248 mm

**B. Randall** CE

NUMBER	SIZE	LENGTH
64-620	3.9 mm	9.9"/252 mm



**Kevorkian Collector™**

NUMBER	LENGTH
66-441	12"/305 mm

**London™**

NUMBER	LENGTH
64-490	11.7"/297 mm

**Endocrette™** CE

NUMBER	LENGTH
66-447	12.3"/312 mm

**C. Sims**

NUMBER	SIZE	LENGTH	CE
392-110	6 mm	11"/280 mm	CE
392-111	7 mm	11"/280 mm	CE
392-112	8 mm	11"/280 mm	CE
392-113	9 mm	11"/280 mm	CE
392-114	11.5 mm	11"/280 mm	CE
392-116	17 mm	11"/280 mm	

## Euro-Med™ Stainless Steel Hooks



### Angle Hooks

NUMBER	DESCRIPTION	LENGTH
36-630	Straight Hook	9.8"/249 mm
36-632	Right Hook	9.6"/244 mm



### Iris Hook

NUMBER	LENGTH
36-625	10.3"/262 mm



### Marlow Corkscrew Hook <sup>CE</sup>

NUMBER	LENGTH
67-451	9.3"/236 mm



### Double IUD Hook <sup>CE</sup>

NUMBER	LENGTH
392-151	10.6"/269 mm



### Emmett Tenaculum Hook <sup>CE</sup>

NUMBER	LENGTH
66-450	9.8"/249 mm



### IUD Hook <sup>CE</sup>

NUMBER	LENGTH
392-148	10.3"/262 mm

## Euro-Med™ Stainless Steel Cervical Tenacula



### Braun Tenaculum Forceps <sup>CE</sup>

NUMBER	LENGTH
64-420	9.9"/252 mm



### Schroeder Tenaculum <sup>CE</sup>

NUMBER	LENGTH
64-410	9.9"/252 mm



### Duplay Tenaculum Double Curved <sup>CE</sup>

NUMBER	LENGTH
61352	11"/280 mm



### Goldstein Grasp Cervical Stabilizer <sup>CE</sup>

NUMBER	LENGTH
64-6026	9.4"/239 mm

# Euro-Med™ Stainless Steel Dressing and Tissue Forceps



**Allis** CE

NUMBER	DESCRIPTION	LENGTH
61981	4 x 5 Teeth	6"
61983	5 x 6 Teeth	7.5"
61882	5 x 6 Teeth	10"



**Potts-Smith Tissue** CE

NUMBER	DESCRIPTION	LENGTH
61871	1 x 2 Teeth	7"
61872	1 x 2 Teeth	8.2"



**Dressing Serrated** CE

NUMBER	LENGTH
61804	5.6"
61805	6.2"



**Foerster Sponge Forceps** CE

Featuring oval tips with serrations, the Foerster is used to hold sponges during the procedure.

NUMBER	LENGTH
16-182	9.7"



**Potts-Smith Dressing** CE

NUMBER	LENGTH
10-202	8.25"



**Adson Dressing Serrated** CE

NUMBER	DESCRIPTION	LENGTH
61830	Serrated Tips	4.75"



**Adson Tissue 1 x 2 Serrated** CE

NUMBER	DESCRIPTION	LENGTH
61864	1 x 2 Serrated Tips	4.75"



**Tissue 1 x 2** CE

NUMBER	DESCRIPTION	LENGTH
61841	1 x 2 Teeth	5.75"
61843	1 x 2 Teeth	7"
61844	1 x 2 Teeth	8"



**Adson Tissue 1 x 2** CE

NUMBER	DESCRIPTION	LENGTH
61863	1 x 2 Teeth	4.75"



**Adson-Brown Grasping** CE

NUMBER	DESCRIPTION	LENGTH
61866	7 x 7 Teeth	4.75"

## Euro-Med™ Stainless Steel Uterine Forceps



**Bozeman Uterine Forceps** €€

NUMBER	LENGTH
392-660	10.2"/257 mm



**Kelly Uterine Forceps** €€

NUMBER	LENGTH
417-131	12.5"/318 mm



**Foerster Uterine Polyp Forceps** €€

NUMBER	DESCRIPTION	LENGTH
415-140	Straight	10"/254 mm
415-141	Curved	10"/254 mm

## Euro-Med™ Stainless Steel Scissors



**Operating Scissors (Straight)** €€

NUMBER	DESCRIPTION	LENGTH
62002	B/B	5.5"
62008	S/B	5.5"
62014	S/S	5.5"
62004	B/B	6.5"
62005	B/B	7.5"
62010	S/B	6.5"
62011	S/B	7.23"



**Spencer Suture Scissors** €€

NUMBER	DESCRIPTION	LENGTH
62164	Straight	3.5"

## Euro-Med™ Stainless Steel Needle Holders



**Mayo-Hegar Needle Holder (Wider Jaw)** €€

NUMBER	LENGTH
62809	5.5"
62810	6.3"
62815	7.1"
62820	8.1"



**Heaney Needle Holder (Curved)** €€

NUMBER	LENGTH
62825	8.3"



**Crile-Wood Needle Holder (Narrow Jaw)** €€

NUMBER	LENGTH
62840	7.1"



**Metzenbaum Scissors (Straight or Curved)** €€

Straight		Curved	
NUMBER	LENGTH	NUMBER	LENGTH
62116	7"	62131	5.9"
62127	8"	62118	7"
62120	9.1"	62119	8.1"
62122	11.2"	62124	9"
		62126	11"



**Mayo Scissors (Straight or Curved)** €€


Straight		Curved	
NUMBER	LENGTH	NUMBER	LENGTH
62063	5.7"	62067	5.5"
62065	6.8"	62069	6.75"
62066	9"	62070	9"



**Iris Scissors (Straight or Curved)** €€

NUMBER	DESCRIPTION	LENGTH
62140	Straight S/S	4.5"
62143	Curved S/S	4.5"

# Endometrial Samplers



## Pipelle® Endometrial Suction Curette CA\* CE

NUMBER	DESCRIPTION
8200	Pipelle Endometrial Suction Curette (25/Box)



## Pipet Curet™ Endometrial Suction Curette CA\* CE


NUMBER	DESCRIPTION
MX140	3 mm (50/Box) – Greater integrity for stenotic Os
MX145	3 mm Soft Yellow (50/Box) – Pliable for molding catheter
MX150	4 mm (50/Box) – Multiparous patient



## Uterine Explora™ Model I Curette with Sharp Randall-Type Cutting Edge and Vacu-Lok Syringe CA\* CE

Disposable Explora Model incorporates a long nylon curette with Randall-type cutting edge, stylet and Vacu-Lok syringe that adjusts suction from approximately 15- to 20-inch Hg negative pressure (6 cc to 12 cc).

NUMBER	DESCRIPTION
MX120	Uterine Explora with Vacu-Lok Syringe, 3 mm, Model I (25/Box)
MX124	Uterine Explora with Vacu-Lok Syringe, 4 mm, Model I (25/Box)



## Uterine Explora™ Model II Curette with Dual Hole Design and Vacu-Lok Syringe CA\* CE

Uterine Explora Model II offers the same features as the Model I plus it has a dual-port opening. The Vacu-Lok syringe allows you to adjust the suction from approximately 15- to 20-inch Hg negative pressure (6 cc to 12 cc).


NUMBER	DESCRIPTION
MX130	Uterine Explora II with Vacu-Lok Syringe Model II (25/Box)



# Dilators

## Comfort-Flex Dilator® Set CA\*

- Os Locator (1 ring) – opens the external Os
- Canal Dilator (2 rings) – traces the natural path of the cervical canal while steadily dilating
- Fundus Sound™ (3 rings) – opens the internal Os and sounds the fundus to determine uterine depth



NUMBER	DESCRIPTION
64-376	Comfort-Flex Dilator Set, 3 dilators per set (Os Locator, Canal Dilator, Fundus Sound)



## Milex® Vaginal-Hymenal Silicone Dilators CA\*

NUMBER	DESCRIPTION
MX20-CA	Milex Vaginal Dilator – Set of 4: Small, Medium, Large, Extra-Large without Trimo-San
MX20I	Milex Vaginal Dilator Single – Small 5 1/8" x 1 5/16"
MX20IM	Milex Vaginal Dilator Single – Medium 5 1/2" x 1 1/16"
MX20IL	Milex Vaginal Dilator Single – Large 5 3/4" x 1 1/4"
MX20IXL	Milex Vaginal Dilator Single – Extra-Large 6 1/8" x 1 1/2"

## Milex® Vaginal-Rectal Silicone Dilators CA\*

NUMBER	DESCRIPTION
MX25-CA	Milex Rectal Dilator – Set of 4: Small Rectal, Medium Rectal, Large Rectal, Extra-Large Vaginal-Rectal without Trimo-San (International sales only)
MX25IS	Milex Vaginal-Rectal Dilator Single – Small 3" x 1 1/16"
MX25IM	Milex Rectal Dilator Single – Medium 3 1/8" x 1 5/16"
MX25IL	Milex Rectal Dilator Single – Large 3 9/16" x 1 1/8"
MX25IXL	Milex Vaginal-Rectal Dilator Single – Extra-Large 3 15/16" x 1 5/16"

## Milex® Vaginal Silicone Dilators CA\*

NUMBER	DESCRIPTION
MX23-CA	Milex Dilator – Set of 4: Small, Medium, Large, Extra-Large without Trimo-San
MX23I	Milex Dilator Single – Small 2 1/16" x 1 1/2" x 3/8"
MX23IM	Milex Dilator Single – Medium 2 5/16" x 9/16" x 7/16"
MX23IL	Milex Dilator Single – Large 2 11/16" x 5/8" x 1 1/2"
MX23IXL	Milex Dilator Single – Extra-Large 3" x 1 1/16" x 9/16"

# Disposable Uterine Sound and Dilators




## Fundus Sound™ XM CA\*

NUMBER	DESCRIPTION
1196	Fundus Sound XM Disposable Uterine Sound (25/Box)



## Milex® Dilator / Sound Set CA\*

NUMBER	DESCRIPTION
MX11	Milex Os Locator (25/Box)
MX12	Milex Canal Dilator (25/Box)
MX13	Milex Fundal Sound (25/Box)
MX21	Milex Dilator/Sound Set, 3 dilators per set: MX11, MX12, MX13 (10/Box)



## Os Finder™ CA\* CE

NUMBER	DESCRIPTION
1176	Os Finder (50/Box)



## Flexisound™ CA\* Intrauterine Sound

Graduated in centimeters.

NUMBER	DESCRIPTION
920008	Flexisound Disposable Intrauterine Sounds (25 per Box)



## Disposable Os Finders

- Disposable Os Finders
- Os Finder
  - Canal Finder
  - Fundus Finder

NUMBER	DESCRIPTION
920019	Disposable Os Finders – Set of 3 (10 Sets per Box)

# Cervical and Endocervical Samplers

# Diagnostic Disposables



## Medscand® Cytobrush Plus® Kits

NUMBER	DESCRIPTION
10115-10/50	Cytobrush Plus Kits – Cell Collector (sets include brush, scraper, fixative, slide and slide transporter); 10 Trays with 50 Pap Smear test sets (500 total Tests/Case)
60115-10/50	Cytobrush Plus Cell Collector – Gentle Touch, 10 Trays with 50 Pap Smear test sets (500 total Tests/Case)
60115-20/50	Cytobrush Plus Cell Collector – Gentle Touch, 20 Trays with 50 Pap Smear test sets (1,000 total Tests/Case)



## Medscand® Cytobrush Plus® CA\* CE

NUMBER	DESCRIPTION
C0021	Cytobrush Plus – 2 Bags, 500 Brushes each (1,000/Box)
C0004	Cytobrush Plus – 10 Bags, 10 Brushes each (100/Box)
C0005	Cytobrush Plus – 25 Bags, 100 Brushes each (2,500/Box)
C0012	Cytobrush Plus Sterile, 1 Brush per Pouch (40/Box)



## Medscand® Cytobrush Plus® GT CA\* CE

NUMBER	DESCRIPTION
C0121	Cytobrush Plus GT – 2 Bags, 500 Brushes each (1,000/Box)
C0104	Cytobrush Plus GT – 10 Bags, 10 Brushes each (100/Box)
C0105	Cytobrush Plus GT – 25 Bags, 100 Brushes each (2,500/Box)
C0112	Cytobrush Plus GT Sterile, 1 Brush per Pouch (40/Box)
C0305	Cytobrush Plus GT Scored – 25 Bags, 100 Brushes each (2,500/Box)



## Medscand® Pap Perfect® Spatula

NUMBER	DESCRIPTION
11080	Pap-Perfect Plastic Spatulas (500/Box) CA* CE



## Milex® Cervical Spatulas

NUMBER	DESCRIPTION
MX600	Milex Cervical Spatulas, Small, Yellow (250/Box)
MX610	Milex Cervical Spatulas, Medium, White (250/Box)



## Papette® Cervical Cell Collector

NUMBER	DESCRIPTION
908001	Papette Cervical Cell Sampler Brush (100 per Bag) CA* CE
908003	Papette Cervical Cell Sampler Brush (10 per Box; 1,000 pieces) CA* CE
908004	Papette Cervical Cell Sampler Brush (40 per Box; 4,000 pieces) CA* CE

## Endocell® Disposable Endometrial Cell Sampler

NUMBER	DESCRIPTION
908014 A	Endocell – 35 per box CA* CE
908015	Endocell – 6 boxes of 35 (210 pieces) CA* CE
908016	Endocell – 15 boxes of 35 (525 pieces) CA* CE



## Word Catheter CA\*


NUMBER	DESCRIPTION
53404	Word Bartholin Gland Catheter. Includes syringe and scalpel. Latex-free. (6/Box)



## Breast Aspiration Dual-Port Biopsy Needle

NUMBER	DESCRIPTION
MX41	Dual-Port Needle with 3-finger control with 12 cc Syringe, 1½" (12/Box) CA*
MX42	Dual-Port Needle, No Syringe, 1½" (5/Box) CA*

# Leisegang™ LM-900 Cryosurgery System



50101



50102



50103



50104



50105



50106



50107



50108



50109

### LEISEGANG LM-900 Probes CA\* CE

NUMBER	DESCRIPTION
50101	Microderm - 2 mm D x 25 mm L
50102	Soll Lash Epilating - 2 mm D x 27 mm L
50103	Condyloma - 5 mm D x 25 mm L
50104	Derm Bevel - 5 mm D x 24 mm L
50105	General Purpose - 19 mm D x 8 mm L
50106	Large Derm Bevel
50107	Curved Anorectal - 8 mm D x 38 mm L
50108	General Purpose - 19 mm D x 8 mm L
50109	Exocervical - 19 mm D x 5 mm L

NUMBER	DESCRIPTION
50501	LM-900 N <sub>2</sub> O Cryosurgery Unit with Case CE
50502	LM-900 CO <sub>2</sub> Cryosurgery Unit with Case CE

### Accessories

NUMBER	DESCRIPTION
50000	Empty 20 lb. N <sub>2</sub> O
50010	Empty 20 lb. CO <sub>2</sub>
50003	20 lb. Cylinder Cart
50008	E-Cylinder Adapter CO <sub>2</sub> CE
50009	Carrying Case
50011	O-Rings (10 μ) CE
50012	Cleaning Plugs (1 μ) CE

# Cryosurgical Systems:

One- and Two-Trigger Systems Provide Conductive, Reusable Cryotips and Disposable Plastic Shields Available for Cryosurgical Treatment



## LL100™ System (N<sub>2</sub>O) Two-trigger

Freeze and defrost only when triggered. Use with cryotips and disposable plastic shields.

## LL100™ System (N<sub>2</sub>O) & LLCO<sub>2</sub>™ System

NUMBER	DESCRIPTION	
900001	LL100 Cryosurgical System – N <sub>2</sub> O (Freezer Only)	
900160	LLCO <sub>2</sub> Cryosurgical System – CO <sub>2</sub> (Freezer Only)	
901060	N <sub>2</sub> O Cylinder (unfilled)	CA* CE
901061	CO <sub>2</sub> Cylinder (unfilled)	CE
901011	Cart for 20 lb. Cylinder	

## LLCO<sub>2</sub>™ System Single-trigger for CO<sub>2</sub> Gas

- Squeeze trigger to freeze; release to defrost
- Use with medical-grade CO<sub>2</sub> gas for optimal performance
- Use with CooperSurgical's standard cryotips with disposable plastic shields

Both the LL100 and the LLCO<sub>2</sub> are available in three-tip package specials for additional savings. Choose three cryotips from our comprehensive offerings when ordering. The following are recommended for general/family practice procedures: T-1905, T-2514, and T-0500 for GYN procedures and the T-0524, T-0219, and T-1900.



## WA1000B™ Cryo Console System

Uses LL100™ two-trigger or LLCO<sub>2</sub>™ single-trigger freezers

- Temperature Indicator with integrated timer
- Cart carries a single 20 lb. cylinder
- Use with CooperSurgical's large selection of cryotips and shields

NUMBER	DESCRIPTION	
900506-1-N2O	WA1000B Cryosurgical Console – Blank N <sub>2</sub> O LL100	CA* CE
900506-1-CO2	WA1000B Cryosurgical Console – Blank CO <sub>2</sub> LL100	CA*
900506-2-N2O	WA1000B Cryosurgical Console – Digital Temp. (N <sub>2</sub> O) Indicator, Timer, Base, and LL100	CE
900506-2-CO2	WA1000B Cryosurgical Console – Digital Temp. (CO <sub>2</sub> ) Indicator, Timer, Base, and LLCO <sub>2</sub>	



## Cryosurgical Tips for Use with LL100™ and LLCO<sub>2</sub>™ CA\* CE

Conductive, thermal, reusable, sterilizable cryotips with disposable plastic shields.

NUMBER	DESCRIPTION	NUMBER	DESCRIPTION	NUMBER	DESCRIPTION	NUMBER	DESCRIPTION
900300AA	(T-0200)	900200AA	(T-0219)	900206AA	(T-0826)	900212AA	(T-1920)
900301AA	(T-0500)	900201AA	(T-0219)	900207AA	(T-0823)	900213AA	(T-2500)
900302AA	(T-0700)	900202AA	(T-0524)	900208AA	(T-1900)	900214AA	(T-2507)
900303AA	(T-0900)	900203AA	(T-0519)	900209AA	(T-1905)	900215AA	(T-2514)
900305AA	(T-1300)	900204AA	(T-0832)	900210AA	(T-1920)	900217AA	(T-0832)
900306AA	(T-1500)	900205AA	(T-0819)	900211AA	(T-1910)	900221AA	Disposable Plastic Shields for Cryotips (100 per Bag)

# H/S Elliptosphere® Catheter Set and Tray

Sonohysterography and Hysterosalpingography

Available in two versions – with and without an integrated stylet – the H/S Elliptosphere features latex-free construction.



61-4205 H/S Elliptosphere Tray



NUMBER	DESCRIPTION	CA*	CE
61-4005	H/S Elliptosphere Catheter, 5 Fr (10/Box)	CA*	CE
61-4005S	H/S Elliptosphere Catheter Set with Integrated Stylet (10/Box)	CA*	CE
61-4205	H/S Elliptosphere Tray with 5 Fr Catheter (10/Case)		

# H/S Catheters

For Visualization of Uterine Pathology

## H/S Catheter Set CA\* CE

Designed for use in sonohysterography and hysterosalpingography, the catheter can be used with oil- or saline-based media depending on the procedure requirements.

- Available with integrated stylet
- Comes in both 5 and 7 French

NUMBER	DESCRIPTION
61-5005	H/S Catheter, Kraton, 5 Fr (10/Box)
61-5007	H/S Catheter, Kraton, 7 Fr (10/Box)
61-0020	H/S Stylet Catheter Stiffener (10/Box)

## H/S Procedure Tray

- Fenestrated drape
- Disposable speculum
- Disposable cervical dilator
- Three swab sticks
- Sterile wrap
- Inner tray for povidone-iodine solution
- One pack of povidone-iodine solution
- 36-in. extension tube
- 20 cc syringe with 18-gauge needle
- Lubricating jelly packet
- Sanitary napkin

NUMBER	DESCRIPTION
61-5205	H/S Procedure Tray – Round 5 Fr
61-5207	H/S Procedure Tray – Round 7 Fr
61-3205	H/S Procedure Tray – Round 5 Fr with Stylet
61-3207	H/S Procedure Tray – Round 7 Fr with Stylet



## HSG Procedure Tray

The HSG Procedure Tray is for use in hysterosalpingography and sonohysterography procedures.


NUMBER	DESCRIPTION
6050T	Procedure Tray

# Fertility Assessment



## SelectMucus® Endocervical Aspirator

NUMBER	DESCRIPTION	CA*
SM301	SelectMucus Endocervical Sampler (25/Box)	CA*



## Tampa Infusion Catheter CA\* CE

For sonohysterography – ultrasound examination of the uterus.

NUMBER	DESCRIPTION
61-2005	Tampa Catheter (20/Box)



 **Paragard**<sup>®</sup>  
T380A  
intrauterine copper contraceptive

by CooperSurgical<sup>®</sup>

## Ordering Paragard

Two simple ways to order the Paragard IUS for your patients:



Paragard Direct<sup>SM</sup> is our “buy and bill” option for ordering Paragard that allows you to take advantage of volume discount pricing.



Create an account



Order online, by fax or by phone



Ship to one or more locations

Visit [ParagardDirect.com](http://ParagardDirect.com) or call 1-877-PARAGARD to learn more.



With the Paragard Specialty Pharmacy<sup>SM</sup> program, your patient’s prescription for Paragard is fulfilled through one of our specialty pharmacies and shipped directly to your office.



Process the claim with your patient’s health insurer



Collect any applicable co-pay for the unit, if needed



Ship Paragard directly to your office

To order through a specialty pharmacy:

1. Visit [ParagardAccessCenter.com](http://ParagardAccessCenter.com) to download our specialty pharmacy request form
2. Complete the form and fax it back to the selected specialty pharmacy



Designed for the collection of cord blood and cord tissue stem cells.

**Ordering a CBR® Newborn Stem Cell Preservation Collection Kit**

Have your patient visit [cordblood.com](http://cordblood.com) or call CBR® (Cord Blood Registry<sup>SM</sup>) directly at 1.888.267.3256 or you can email [providersupport@cordblood.com](mailto:providersupport@cordblood.com) to request collection kits for your office or hospital.

# Fetal Monitor



**F9 Fetal Monitor**

Fetal monitor with 2 ultrasound transducers, 1 toco transducer, 1 event marker, 2 reusable belts, ultrasound gel, printer paper, power cord, system user manual.

Optional: CooperSurgical Insight™ software, RS232 and ethernet cables, rechargeable lithium-ion battery.

NUMBER	DESCRIPTION
MS9-108022-UL-R	CooperSurgical Fetal Monitor F9
MS9-108022-SWP	CooperSurgical Insight Software Package (Includes software and cables)
M21-64072	Rechargeable Lithium-ion Battery
MS9-107634	Deluxe Cart
MS9-108166	Base Cart with Basket and Laptop Tray
MS9-108108	Base Cart with Basket
MS3-31528-A1	Transducer 1 (pink, 4-pin)
MS3-31528-A1-R	Transducer 1 (pink, 6-pin)
MS3-107705-A1	Transducer 2 (yellow, 4-pin)
MS3-107705-A1-R	Transducer 2 (yellow, 6-pin)
MS3-31527-A1	TOCO Probe
MS3-31112	Event Marker
MS1R-02264	Transducer Belt (reusable)
M25R-75111	Recording Paper (90 mm x152 mm x 150 p, US standard) (Pack of 3)



**Fetal2EMR™**

Fetal Monitor with the capability of transferring NSTs to electronic medical records.

NUMBER	DESCRIPTION
902300-R	Fetal2EMR Monitor and Accessories
902320	Fetal2EMR Insight™ Software; Includes Installation CD and Connecting Cables
902302	FHR/Ultrasound Transducer (4-pin)
902302-R	FHR/Ultrasound Transducer (6-pin)
902304	Thermo Paper 90 mm x 152 mm x 150 p (Pack of 3)
902340	Fetal2EMR Rolling Cart with Tray
902341	Fetal2EMR Rolling Cart without Tray

# Obstetrical Dopplers and Probes



### LifeDop® 250 Series

Handheld Dopplers with large displays for indication of heart rate and battery status. Rechargeable batteries and audio recording options are available on specific models. Compatible with all LifeDop non-directional probes.

NUMBER	DESCRIPTION	
L250-SD2	LifeDop 250 Disp. with 2 MHz Probe	CA* CE
L250-SD3	LifeDop 250 Disp. with 3 MHz Probe	CA* CE
L250-SDW	LifeDop 250 Disp. with 2 MHz Waterproof Probe	CA* CE
L250R-SD2	LifeDop 250 Disp. with Recharge & 2 MHz Probe	CA* CE
L250R-SD3	LifeDop 250 Disp. with Recharge & 3 MHz Probe	CA* CE
K200	Roll Stand with Storage Basket for L150 or L250	



### LifeDop® 350 Series

The LifeDop 350 Doppler is a tabletop Doppler system designed to detect the fetal heart rate and monitor blood flow. It has a large heart rate display, comes with rechargeable batteries, and accepts all LifeDop non-directional probes.

- Two probe holders for multiple applications
- User-accessible batteries
- Stand and wall mounting available

NUMBER	DESCRIPTION	
L350R-SD2	LifeDop Tabletop Doppler with 2 MHz Probe	CA* CE
K310	Wall Mount for L350R	
K220	Roll Stand with Storage Basket for L350R	



### LifeDop® Non-Directional Obstetrical Probes

All LifeDop non-directional probes are completely interchangeable for use on the L150 and the L250 and L350 series of Dopplers. Most models have user-replaceable cables.

**2 MHz OB Probe** Wide-angle beam for deep penetration needed in late term pregnancy and larger patients.

**3 MHz OB Probe** The general-purpose obstetrical probe for fetal heart detection, from earliest stages through delivery.

**2 MHz Waterproof OB Probe** Waterproof design for water-assisted labor and delivery applications.

NUMBER	DESCRIPTION	
SD2	2 MHz OB Probe	CA* CE
SD3	3 MHz OB Probe	CA* CE
SDW	2 MHz Waterproof OB Probe	CA* CE

# Vascular Assessment Dopplers and Probes



### LifeDop® 150 Series

These handheld, non-display Dopplers for vascular assessment offer full line of non-directional LifeDop probes. Rechargeable batteries and audio recording options are available on certain models.

NUMBER	DESCRIPTION
L150R-SD4	LifeDop 150 with 4 MHz Probe, Rechargeable
L150R-SD5	LifeDop 150 with 5 MHz Probe, Rechargeable
L150R-SD8	LifeDop 150 with 8 MHz Probe, Rechargeable
L150R-STN	LifeDop 150 with 8 MHz Sterilizable Probe, Rechargeable
STNPK1	L150 Doppler with Sterilizable Probe and Stand



### LifeDop® 350

The LifeDop 350 has a large heart rate display, rechargeable batteries, and accepts all LifeDop non-directional probes. Enhanced features include:

- Two probe holders for multiple applications
- User accessible batteries
- Wall mounting

NUMBER	DESCRIPTION
K310	Wall Mount for Added Convenience and Security



### LifeDop® Vascular Assessment Probes

CA\* CE

LifeDop probes are interchangeable for use on the L150 Series and L350 Dopplers, with the exception of the sterilizable probes, which have some restrictions. Most models have user replaceable cables.

- 4 MHz Vascular Probe
- 5 MHz Vascular Probe
- 8 MHz Vascular Probe
- 8 MHz Sterilizable Probe

NUMBER	DESCRIPTION
SD4	4 MHz Vascular Probe
SD5	5 MHz Vascular Probe
SD8	8 MHz Vascular Probe
STN	8 MHz Sterilizable Vascular Probe

# Full-Featured ABI Systems for Peripheral Arterial Disease Diagnosis

## Vista AVS®

### A Full-Featured ABI/Vascular System

The Vista AVS offers three modalities for obtaining systolic pressures and waveforms and is designed to accommodate both standard and custom protocols. The accompanying compact printer prints two ankle waveforms along with the ABI results. The system includes:

- 3 modalities: PVR, bi-directional Doppler probe and PPG probe
- Automatic ABI, TBI and segmental calculation
- Performs seated ABI on mobility impaired patients
- Software for full-page reporting, exam storage and interface with EMR
- 9 blood pressure cuffs: (2) 10 cm, (4) 12 cm, (2) 17 cm contoured thigh cuffs, (1) digit cuff
- Automatic inflation/deflation
- Graphic waveform display with step-by-step ABI instructions

### Optional Features for the Vista AVS®


#### 5 MHz Bi-directional Probe

For deeper vessel penetration on larger patients.

#### Cuff Selector

Up to 10 vascular cuffs can be connected to this unit at the same time, for upper and lower extremity or digit pressure measurements.

NUMBER	DESCRIPTION	
L500VA	Full-featured ABI/Vascular System	CA*
SD5V	5 MHz Bi-directional Probe	
K290	Cuff Selector	
G150	Ultrasound Gel (60 grams)	
VISTCAP1	Extended one year warranty for Vista system	



## Vista ABI™

This system includes:

- 2 modalities: PVR and bi-directional Doppler probe
- Automated cuff inflation/deflation
- Built-in printer
- Graphic waveform display with step-by-step ABI instructions
- 4 cuffs: (2) 10 cm, (2) 12 cm

### Optional Features for the Vista ABI™

#### PPG Package

A sensitive arterial PPG probe for obtaining both toe and limb pressures. Comes with Velcro to secure placement on the digit. This package includes a PPG probe and a 1.9 cm digit cuff for performing TBI exams on patients with non-compressible arteries.

#### PC Download Software

An optional software package for the Vista ABI with full-page reporting and exam storage includes download capability for integration with EMR systems (.pdf and .csv file formats – USB Interface to PC).

NUMBER	DESCRIPTION	
L450VA	Vista ABI System	CA*
L450PG	PPG Package	
L450SF	PC Download Software	
SD5V	5 MHz Bi-directional Probe	
G150	Ultrasound Gel (60 grams)	

\*The Vista ABI does not accommodate segmental exams.



# Portable and Handheld ABI Systems



## LifeDop® 300 ABI System with Four Cuffs and Aneroid

A Handheld ABI System with ABI Calculation and Dual Waveform Printout

Portable bi-directional Doppler with keypad for entering pressure data and calculating the indices. The accompanying compact printer prints two ankle waveforms along with the ABI results.

System includes:

- Printer
- 8 MHz bi-directional probe
- Adhesive-backed label paper
- ABI chart, forms
- ABI video, examination guide, quick reference guide, and ABI glossary
- Carrying case



## LifeDop® 250 ABI System with Cuffs and Aneroid

Portable bi-directional Doppler with compact printer. The system includes handheld Doppler for diagnosing PAD, cuffs and aneroid package, printer to capture waveforms, an 8 MHz bi-directional probe, adhesive-backed label paper, carrying case, and ABI reference documents and video.

NUMBER	DESCRIPTION	
L300AC	LifeDop 300 ABI System with 4 Cuffs and Aneroid	CA* CE
L250AC	LifeDop 250 ABI System with 4 Cuffs and Aneroid	CA* CE
L250AB	LifeDop 250 ABI System with 2 Cuffs and Aneroid	CA* CE
SD5B	5 MHz Bi-directional Probe	
SD8B	8 MHz Bi-directional Probe	
K180	Adhesive-backed label paper	
G150	Ultrasound Gel	



## Vantage ABI™

The Vantage ABI™ is a cuff-based system that uses Digital Fourier Oscillometric or DFO waveform analysis. The accompanying compact printer prints two ankle waveforms along with the ABI results.

### VANABI – Vantage ABI Package 1

- Vantage ABI
- Graphic display with convenient touch-screen to navigate the system
- USB memory stick for data transfer and storage (EMR interface)
- Four blood pressure cuffs: (2) 10 cm cuffs for ankles R/L and (2) 12 cm cuffs for arms R/L
- Automatic inflation/deflation

### VANMAX – Vantage ABI Package 2

- Vantage ABI
- Graphic display with touch-screen to navigate the system
- USB memory stick for data transfer and storage (EMR interface)
- Built-in printer for printing ankle waveforms on adhesive-backed label paper
- Four blood pressure cuffs: (2) 10 cm cuffs for ankles R/L and (2) 12 cm cuffs for arms R/L
- Automatic inflation/deflation

NUMBER	DESCRIPTION
VANABI	Vantage ABI – Package 1
VANMAX	Vantage ABI with Stand – Package 2
K180	Adhesive-backed label paper
C019	Cuff and Hose Assembly – 1.9 cm Digit Cuff (VISTA Systems Only)
C100	Cuff and Hose Assembly – 10 cm
C101	Cuff and Hose Assembly – 10 cm – Ankle Left
C102	Cuff and Hose Assembly – 10 cm – Ankle Right
C120	Cuff and Hose Assembly – 12 cm
C121	Cuff and Hose Assembly – 12 cm – Arm Left
C122	Cuff and Hose Assembly – 12 cm – Arm Right
C170	Cuff and Hose Assembly – 17 cm Contoured Thigh Cuff

# Ophthalmic Cryosurgery

## Insulating Sleeves

NO.	DESCRIPTION
P805	For Probes 124/125/142/161/162/166/167 (Retinal) (6/Pkg)
805-1	For Probes 143/144 (Hammerhead) (6/Pkg)
32709	For Probes 128/132/147 (Micro) (1/Pkg)
32710	For Probes 126/127 (Cataract) (6/Pkg)
32787	For Probes 139/163 (Retinal Spherical) (1/Pkg)
33098	For Probes 130/164/165 (Glaucoma) (6/Pkg)
33555	For Probes 120/121 (Retinal Pediatric) (1/Pkg)



# Wallach® UltraFreeze™

Portable cryosurgical system for most office skin procedures. The system includes an assortment of luer lock open spray apertures. A plastic reusable plate with six conical openings is available to help direct liquid nitrogen spray in concentrated areas. Available in .5 liter.

## Wallach UltraFreeze™

NUMBER	DESCRIPTION
900108	Accu-Shield Cone Spray Plate (set of 6 varied cones)
900111-1	Manual Discharge Device for 10 liter Dewar
900111-3	Manual Discharge Device for 20 liter Dewar
900112	UltraFreeze Stand
900076	UltraFreeze Liquid Nitrogen Sprayer (.5 liter)*
900113	Aperture kit for UltraFreeze (4 spray apertures and one 20 gauge bent needle)



# Index

<b>Biopsy Punches</b>		<b>H/S Catheters</b> .....69	<b>Obstetrical Care</b> .....73-75
Euro-Med.....36-39		<b>H/S Elliptosphere Catheter</b> .....68	<b>Ophthalmic Cryosurgery</b> .....82
<b>Cervical Samplers</b> .....62-63	<b>Euro-Med Instruments</b> <b>(see also LEEP)</b>	<b>Hooks</b>	<b>Pessaries:</b> <b>See Urinary Incontinence</b>
<b>Colposcopy</b>	Biopsy Punches.....36-39	LEEP.....24-25	<b>Paragard</b> .....70-71
OptiK Model 1.....29	Needle Holders.....56	Stainless Steel.....52	<b>Pipelle</b> .....58
OptiK Model 2.....30-31	Sounds and Dilators.....48-49		<b>Pipet Curet</b> .....58
Wallach® Zoom Magnification Colposcopes.....32	Stainless Steel Curettes.....50-51	<b>Hysteroscopy</b>	<b>Retractors</b>
<b>Cord Blood Preservation</b>	Stainless Steel Endocervical Specula.....44-45	Endosee® Advance.....6-7	LEEP.....22
CBR.....72	Stainless Steel Forceps.....54-56	<b>Lateral Wall Retractors</b>	Stainless Steel.....42, 46
<b>Cryosurgery</b>	Stainless Steel Hooks.....52	LEEP.....22	<b>Scissors – Stainless</b> .....57
Leisegang™ LM-900	Wall Retractors.....46	Stainless Steel.....46	<b>Sounds</b> .....48, 61
Cryosurgery System.....64-65	Stainless Steel Scissors.....57	<b>LEEP</b>	<b>Specula</b>
Cryosurgical Systems.....66-67	Stainless Steel Specula.....40-43	Accessories.....18-19	Dimensions.....42-43
Ophthalmic Cryosurgery.....82	Stainless Steel Tenacula.....53	Electrodes.....14-15	LEEP.....20-21
UltraFreze.....83	<b>Fertility and Family Planning</b>	Endocervical Specula.....25	Stainless Steel.....40-41
<b>Cytobrush</b> .....62	Fertility Assessment.....68	Fischer Cone Biopsy Excisor.....16	ER-Spec™.....34
<b>Curettes-Endocervical</b>	H/S Catheter.....69	Forceps.....24	Her-Spec®.....34
Stainless Steel.....50	H/S Elliptosphere Catheter.....68	Hooks.....24-25	OfficeSpec®.....35
<b>Curettes-Endometrial</b>	<b>Fetal Monitors</b>	Lateral Wall Retractors.....22	Anospec®.....35
Disposable.....58	F9.....73	LEEP PRECISION™	<b>Tenacula</b>
Stainless Steel.....51	Fetal 2EMR.....73	Integrated System.....12-13	LEEP.....24
<b>Dilators</b>	<b>Forceps</b>	Procedure Essentials.....17	Stainless Steel.....53
Milex Vaginal-Hymenal	LEEP.....24	Retractors.....22	<b>Tubing</b>
Silicone Dilators.....60	Stainless Steel.....54-56	LEEP.....25	Smoke Evacuator Tubing.....18
Milex Vaginal-Rectal Silicone	<b>General Examination:</b>	Tenacula.....24	Tubing Set.....19
Dilators.....60	<b>Diagnostic Disposables</b>	Tubing.....18-19	<b>Urinary Incontinence</b>
Milex Vaginal Silicone Dilators.....60	Breast Aspiration Dual-Port	Vaginal Specula.....20-21	Milex Pessaries.....8-11
Milex Dilator/Sound Set.....61	Biopsy Needle.....63	Wallach® Quantum 2000	<b>Uterine Infusion Catheters</b>
<b>Dopplers</b>	Cervical and	Accessories.....26-27	H/S Catheters.....69
Obstetrical Doppler Probes	Endocervical Samplers.....62-63	<b>LEISEGANG</b>	H/S Elliptosphere Catheter.....68
and Accessories.....74-75	Milex Cervical Spatulas.....62	Colposcopy.....28-31	<b>Uterine Sounds</b>
Vascular Dopplers.....76-81	Word Catheter.....63	LM-900 Cryosurgery System.....64-65	Stainless Steel.....48
<b>Electrodes: See LEEP</b>	<b>Gynecological Specialty</b>	<b>Lighted Specula – Vaginal</b>	Disposable.....61
<b>Endocervical Samplers</b>	<b>Instruments</b>	and Anal	
Cytobrush.....62	Euro-Med Stainless Steel	ER-Spec™.....34	
<b>Endocervical Specula</b>	Dilators.....48-49	Her-Spec®.....34	
LEEP.....25	Endometrial Curettes.....51	OfficeSpec®.....35	
Stainless Steel.....44-45	Euro-Med Stainless Steel	Anospec®.....35	
<b>Endometrial Samplers</b>	Endocervical Curettes.....50	<b>Milex</b>	
Disposable.....63	Euro-Med Stainless Steel	Dilators.....62	
EMB Samplers.....58	Endocervical Specula.....44-45	Pessaries.....8-11	
<b>Endosee® Advance</b> .....6-7	Euro-Med Stainless Steel Hooks.....52	<b>Needle Holders</b> .....56	
	Euro-Med Stainless Steel	<b>Needles</b>	
	Specula.....40-43	Potocky Needle.....17	
	Forceps.....54-56	Biopsy Needle.....63	
	Euro-Med Lateral Wall Retractors.....46		
	Euro-Med Stainless Steel		
	Uterine Sounds.....48		

## Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers' expectations, but industry standards as well. CooperSurgical's Customer Service Department provides full-service order entry with extended hours to facilitate order placement. As an added convenience we also now offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

☑ Available for sale where CE mark is required

CA Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most up-to-date information.

## Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products and consumables, will be free from defects in materials and workmanship for one year from the date of purchase.

CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.


THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.


COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

## 5 Ways to Order

 **Local Sales Representative**  
www.coopersurgical.com/pages/salesrepsearch.aspx

 **Email**  
U.S.: Orders@coopersurgical.com  
Int'l: InternationalCustomerservice@coopersurgical.com

 **E-Ordering (U.S. Only)**  
Call 800.243.2974 to set up your account today.

 **Phone Orders**  
800.243.2974 (8:00 AM – 8:30 PM EST)  
Int'l: 203.601.9818

 **Fax Orders**  
800.262.0105 (24 hours a day)  
Int'l: 203.601.4747

## Returned Goods Policy

- Goods will be accepted for return for the following reasons:**
  - If shipment was made without Customer's authorization or order
  - If incorrect goods were shipped
  - If defective goods are covered by the standard warranty
- To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without an RMA number.  
**Please have the following information:**
  - Reason you wish to return the goods
  - Quantity, description, part number, serial number of the goods
  - Date of Receipt of Order
  - Customer's Purchase Order and the CooperSurgical Invoice Number**All used instruments must be cleaned and sterilized prior to being shipped for return.**
- Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.
- If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:
  - There is no restocking fee if goods are returned within 30 days of shipment
  - There is a restocking fee of 25% if goods are returned before 31-60 days of shipment
  - Goods must be returned unused, in the original carton, and in marketable condition
  - Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
  - No return after 60 days
  - No refund on sterile, single-use disposable products

Send All Returns to: **CooperSurgical**  
**95 Corporate Drive**  
**Trumbull, CT 06611**

## Trademark Notices

All trademarks are the property of their respective owners.

Endosee® Advance

Milex®

LEEP PRECISION™

Fischer Cone Biopsy Excisor®

Pipelle®

Leisegang™

ER-Spec®

Her-Spec®

Euro-Med®

H/S Elliptosphere®

Paragard® Intrauterine  
Copper Contraceptive

Cord Blood Registry<sup>SM</sup>

Wallach®



CooperSurgical®

800.243.2974

[www.coopersurgical.com](http://www.coopersurgical.com)

95 Corporate Drive, Trumbull, CT 06611